

# High risk PCI & Hemodynamic support (IABP, Impella, ECMO)

**Adam Taleb, MD FACC**

Interventional Cardiologist

Athens Medical Center



**ΟΜΙΛΟΣ ΙΑΤΡΙΚΟΥ  
ΑΘΗΝΩΝ**

Ιατρικό Κέντρο Αθηνών

# Conflict of interest:


Proctoring for IMPELLA ABIOMED



YOU SHALL NOT PASS!


# High Risk PCI - CHIP

# Definition of the CHIP Population: Complex Higher-Risk (and Indicated) Patients



**Are We Treating These Patients Effectively?**

# Severely Calcified Lesions are Undertreated



12%  
Severely  
Calcified  
Lesions

<5% Total  
Atherectomy\*  
Penetration


\* Total Atherectomy includes: Diamondback 360® Orbital Atherectomy, Rotablator™ and Laser Atherectomy systems

12% severe calcium  
Bourantas CV, et al. Heart. 2014;100:1158-1164. 20% severe Ca++ in DES trials  
Généreux P, et al. J Am Coll Cardiol. 2014;63:1845-54. 6% severe Ca++ in ACS  
Feldman DN, et al. Circulation. 2013;127:2295-2306. 60% of PCIs ACS, 40% stable → 12% severe Ca++

Diamondback 360 is a registered trademark of Cardiovascular Systems, Inc (CSI).  
Rotablator is a filed for trademark of Boston Scientific.

IABP


# Φυσιολογική λειτουργία IABP



# DATA- IABP-SHOCK

## Shock and Acute MI


- Randomized, prospective, trial
- 600 patients with NSTEMI/STEMI
- Shock (SBP <90mmHg for >30min, pulmonary congestion, end-organ hypoperfusion)
- Early PCI expected
- PEP: 30d mortality
- IABP vs no IABP



- 30 crossovers to IABP group (26 protocol deviations, not done for mechanical complications)
- **Selection bias:** Rapidly deteriorating patients may not have been enrolled; biasing the study population towards patients exhibiting signs of stability with vasopressor and inotropic support.



# Counterpulsation Requires Native Left Ventricular Pulsation



The more dysfunction the ventricle,  
the *less functional* an IABP becomes

# Recommendations for the use of short-term mechanical circulatory support in patients with cardiogenic shock

| Recommendations  | Class      | Level    |
|--|------------|----------|
| Short-term MCS should be considered in patients with cardiogenic shock as a BTR, BTD, BTB. Further indications include treatment of the cause of cardiogenic shock or long-term MCS or transplantation.            | <b>IIa</b> | <b>C</b> |
| IABP may be considered in patients with cardiogenic shock as a BTR, BTD, BTB, including treatment of the cause of cardiogenic shock (i.e. mechanical complication of acute MI) or longterm MCS or transplantation. | <b>IIb</b> | <b>C</b> |
| IABP is not routinely recommended in post-MI cardiogenic shock.  | <b>III</b> | <b>B</b> |

BTB = bridge to bridge; BTD = bridge to decision; BTR = bridge to recovery; IABP = intra-aortic balloon pump; MCS = mechanical circulatory support; MI = myocardial infarction.

# CASE

- 51 M with PMH of ischemic cardiomyopathy LVEF 20%, CAD s/p PCI in LAD
- Presented after ICD shock x1 for VF
- Cath in OSH

Im: 1/53

Se: 6

1/1/1971 M  
G.H.A Alexandra  
R202307311001383  
Cardiac  
Left Coronary 15 fps



Im: 1/53

Se: 1

1/1/1971 M  
G.H.A Alexandra  
R202307311001383  
Cardiac  
Left Coronary 15 fps



Im: 1/56

Se: 3

1/1/1971 M  
G.H.A Alexandra  
R202307311001383  
Cardiac  
Left Coronary 15 fps

Im: 1/57

Se: 4

1/1/1971 M  
G.H.A Alexandra  
R202307311001383  
Cardiac  
Left Coronary 15 fps

# Decision for IABP

Im: 1/52  
Se: 1

17/4/1971 M  
IATRIKO KENTRO ATHINON  
6945222.7483  
CORONARIES  
CORONARIES

Im: 1/5  
Se: 1

17/4/1971 M  
IATRIKO KENTRO ATHINON  
6945222.7483  
CORONARIES  
CORONARIES




Im: 1/6  
Se: 1

17/4/1971 M  
IATRIKO KENTRO ATHINON  
6945222.7483  
CORONARIES  
CORONARIES



17/4/1971 M  
IATRIKO KENTRO ATHINON  
6945222.7483  
CORONARIES  
CORONARIES







## 2<sup>nd</sup> year follow up

- Pt is tolerating maximum heart failure medical therapy
- He is NOT on the transplant list any more, since LVEF recovered to 35 – 40 %

IMPELLA και iVAC

# iVAC 2L





# IMPELLA

- Continuous axial-flow pump
- Non-pulsatile, no ECG or waveform triggering
- From LV to aorta
- IMPELLA CP, 5.0, 5.5
  
- From IVC to PA
- IMPELLA RP



# Axial Pump Flow Principles

Pump flow is a function of:

– The speed of the impeller

↑ Speed → ↑ Flow

↓ Speed → ↓ Flow

AND


– The difference in pressure across the pump

↑ Pressure gradient → ↓ Flow

↓ Pressure gradient → ↑ Flow

At any given speed, increased pressure gradient across the aortic valve will decrease flow

# The Spectrum of Patients Requiring Acute MCS Therapy



Condition 1: 'Normal'

Condition 2: AMI

Condition 3: Compensated HFrEF

Condition 4: Cardiogenic Shock (AMI or HFrEF)




ORIGINAL ARTICLE

### Cardiogenic Shock From Heart Failure Versus Acute Myocardial Infarction: Clinical Characteristics, Hospital Course, and 1-Year Outcomes


Shashank S. Sinha, MD, MSc<sup>\*</sup>, Carolyn M. Rosner, MSN, NP-C, MBA<sup>†</sup>, Behnam N. Tehrani, MD, Aneel Maini, BS, Alexander G. Truesdell, MD<sup>‡</sup>, Seiyon Ben Lee, PhD, Pramita Bagchi, PhD, James Cameron, MS, Abdulla A. Damluji, MD, PhD, MPH<sup>§</sup>, Mehul Desai, MD, Shashank S. Desai, MD, MBA<sup>¶</sup>, Kelly C. Epps, MD, Christopher deFilippi, MD<sup>||</sup>, M. Casey Flanagan, MD<sup>||</sup>, Leonard Genovese, MD<sup>||</sup>, Hala Moukhachen, MD, James J. Park, BA, Mitchell A. Psotka, MD, PhD, Anika Raja, BS, MPH, Palak Shah, MD, MS<sup>||</sup>, Matthew W. Sherwood, MD, MHS<sup>||</sup>, Ramesh Singh, MD<sup>||</sup>, Daniel Tang, MD, Karl D. Young, PA<sup>||</sup>, Timothy Welch, MD, Christopher M. O'Connor, MD, and Wayne B. Batchelor, MD, MHS<sup>||</sup>


“the acute onset seen in the AMI-CS phenotype may be associated with less tolerance of the hemometabolic perturbation. In contrast, patients with HF-CS appear to have an acute on chronic presentation with a degree of chronic compensation”




# Case


- 72 man referred for PCI to distal RCA and LM ISR CTO
- PMH: CAD s/p STEMI 2015 with LM PCI, ischemic HFrEF 20%, hepatitis C with cirrhosis c/b HCC s/p TACE, stroke, COPD
- TTE: LVEF 20%, anterior septal/apical and distal inferior HK, mild AI
- Cardiac MRI showed viable anterior wall
- Turned down by CT surgery as poor surgical candidate
- Meds: aspirin, ticagrelor, crestor, metoprolol succ, Aldactone












- RHC
- IMPELLA placement
- PCI of the RCA











# IVUS






After stenting  
the mid LAD...











Deploy the  
LM/LAD stent









# Ostial flash balloon






ECMO



# ECMO



VA ECMO is composed of:

- Inflow cannula
- Tubing
- Pump
- Oxygenator
- Heat-exchanger
- Outflow cannula









VA ECMO

Heart




# Hemodynamic Effects of MCS Therapy




# THE NEED FOR ACTIVE UNLOADING IN CS

- Further Impairment of Vulnerable LV
- Retrograde Perfusion
- LV Thrombus Formation



# V-A ECMO LOADS THE LEFT VENTRICLE



\* Under study

Fincke R, et al. J Am Coll Cardiol 2004  
 den Uil CA, et al. Eur Heart J 2010  
 Mendoza DD, et al. Am Heart J 2007  
 Torgersen C, et al. Crit Care 2009  
 Torre-Amione G, et al. J Card Fail 2009


Suga H. Am J Physiol 1979  
 Suga H, et al. Am J Physiol 1981  
 Burkhoff D, et al. Am J Physiol Heart Circ Physiol 2005  
 Burkhoff D. Mechanical Properties Of The Heart And Its Interaction With The Vascular System. (White Paper) 2011

Sauren LDC, et al. Artif Organs 2007  
 Meyns B, et al. J Am Coll Cardiol 2003  
 Rimmelink M, et al. Catheter Cardiovasc Interv 2007  
 Aql RA, et al. J Nucl Cardiol 2009  
 Lam K, et al. Clin Res Cardiol 2009

Reesink KD, et al. Chest 2004  
 Esposito M, et al. J Am Coll Cardiol 2018  
 Rimmelink M, et al. Catheter Cardiovasc Interv 2010  
 Naidu SS. Circulation 2011  
 Weber DM, et al. Cardiac Interventions Today Supplement Aug/Sep 2009



# EFFECTS OF LV DISTENTION DURING V-A ECMO

Myocardial Recovery  
Inversely Related to the Degree of LV Distention (LVD)  
N=121




VA ECMO

Heart



# MCS complications



**Fig. 1** Complications of MCS by category. Each *color* represents a unique complication category. Complications across more than one category are shown in a combination of colors of the two categories







ΟΜΙΛΟΣ ΙΑΤΡΙΚΟΥ  
ΑΘΗΝΩΝ

Ιατρικό Κέντρο Αθηνών



: [taleb.adam@gmail.com](mailto:taleb.adam@gmail.com)



: @adam\_ta1eb

*Ευχαριστώ για την προσοχή σας*